

Ecoship Box

ALL SINGLE USE



1

Open Ecoship Box, ensure paperwork and empty box are stored in a safe place. Distribute sharp containers throughout facility until full.



2

Once containers are full, ensure all sharp containers are locked.



3

Ensure the inner flaps are upright when placing locked sharp containers inside the box. This will ensure the box stays secure through transport.



4

Place locked sharp containers in standing position inside the Ecoship Box.



5

To close the box, fold side flaps over closed sharp containers.



5a

Fold lid and tuck into the front of box. **Ensure security tab is fully inserted into lid.**



6

If applicable, place return bill of lading inside the provided labelope. If you were not provided a return bill of lading, contact Daniels to have one sent to you via email or fax.



7

Contact Daniels by emailing your service request form to your local Daniels service location. A Daniels or Purolator courier will be sent out to collect your waste and provide new boxes.