

NEW AND IMPROVED

D4 SecureSmart Stainless Bracket



MAKING HEALTHCARE SAFER
SOINS DE SANTÉ PLUS SÉCURITAIRES



Keep your sharps containers secure from tampering.

Materials

- Coverall door is type 304 18 ga. stainless steel all welded construction. Rear mounting plate is extra heavy duty 14 ga. powder coated steel.

Features

- Easy loading with front pivot hinged door (Stainless steel rivets).
- Premium tumbler keyed lock (reversible key).
- Radius and bevelled corners provide a sleek user friendly look.
- Container is well recessed with a tube inlet to prevent unauthorized access.
- Accepts Daniels Health D4 sharps container.

Dimensions: 12" (H) x 5.5" (W) x 5.5" (D)

Product Code: 7124

**Partner it with the
Ecoship D4ES Kit**

