



Pharmasmart

Designed with in-built tamper-proof locks which, when engaged, can only be released at a Daniels processing facility; the Daniels Pharmasmart range of collectors. eliminates the common security risks associated with most pharmaceutical waste containers.



WIDE OPENING

One-hand, one-step deposition of sharps. The optimal access geometry of the lid opening, designed with a wide aperture for large and awkward sharps, minimises disposal-related injuries and restricts hand access to the contained waste.




NO LEAKS

A fixed liquid seal around the rim of the container lid eliminates leakage of bodily fluids during transport. Container leakage is a common hazard with sharps containers that are not equipped with a leakproof seal. Neoprene gasket in lid eliminates opportunity for blood or body fluid exposure during post-closure handling and transport



CANNOT BE REOPENED

Side locks are tamper proof when engaged. This safety feature prevents unauthorised access to the container contents on site and in transport. Temporary closure allows safe and easy transport within healthcare facility.



MAXIMUM STRENGTH HARDENED PLASTIC

The Medical grade hardened plastic casing of the Sharpsmart container is completely impenetrable by the contained sharps. The Sharpsmart container has been test-proven to an amazing 500 cycles proving its durability, sterility and re-usability.



ELIMINATES OVERFILLING

Once the container is full, the in-built safety tray is activated into an upright closed position to restrict access and overfilling. Passive overfill protection prevents staff from overfilling container and eliminates risk of overfilling injuries

Pharmasmart Pharmaceutical Range

REGULAR



P14

14 Litre Regular Pharmasmart designed for collection of small quantities of expired and residual pharmaceuticals that require controlled disposal

Fill Capacity: 6.5 Litres

H 285mm x W 398mm x D 210mm



P22

22 Litre Regular Pharmasmart designed for collection of moderate quantities of expired and residual pharmaceuticals that require controlled disposal

Fill Capacity: 14.5 Litres

H 446mm x W 398mm x D 210mm



P32

32 Litre Regular Pharmasmart designed for collection of larger quantities of expired and residual pharmaceuticals that require controlled disposal

Fill Capacity: 23.5 Litres

H 618mm x W 398mm x D 210mm

ACCESS PLUS



P32 AP

32 Litre Access Plus Pharmasmart provides unrestricted deposit. For use in secure areas such as pharmacies, labs and operating rooms that are secure from public access

Fill Capacity: 23.5 Litres

H 618mm x W 398mm x D 210mm



P64 AP

64 Litre Access Plus Pharmasmart provides unrestricted deposit. For use in secure areas such as pharmacies, labs and operating rooms that are secure from public access

Fill Capacity: 47 Litres

H 726mm x W 398mm x D 354mm